

extendSURE™

CONTROL**LOT SPECIFIC VALUE ASSIGNMENT CARD****LEVEL 1****Haemoglobin A1c****UNITS****MEAN****RANGE**

NGSP ALIGNED

%

6.1

(5.1-7.1)

IFCC ALIGNED

mmol/mol

44

(37-51)

g/dL

0.60

(0.50-0.70)

mmol/L

0.37

(0.31-0.43)

Total Haemoglobin

g/dL

13.97

(11.59-16.35)

mmol/L

8.66

(7.19-10.13)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4197



IVD

EC REP

Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New ZealandEMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

USE BY 2021-12-31

Version 10. Revision date 07/19

extendSURE™

CONTROL**LOT SPECIFIC VALUE ASSIGNMENT CARD****LEVEL 2****Haemoglobin A1c****UNITS****MEAN****RANGE**

NGSP ALIGNED

%

11.7

(9.9-13.5)

IFCC ALIGNED

mmol/mol

104

(88-120)

g/dL

1.46

(1.24-1.68)

mmol/L

0.91

(0.77-1.05)

Total Haemoglobin

g/dL

13.89

(11.47-16.31)

mmol/L

8.61

(7.11-10.11)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

BECKMAN COULTER DxC CHEMISTRY SYSTEMS

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4197



IVD

EC REP

Canterbury Scientific Ltd
71 Whiteleigh Avenue
Addington, 8011
Christchurch
New ZealandEMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

USE BY 2021-12-31

Version 10. Revision date 10/19

Printed to either side of
225g card