

# extendSURE™

**CONTROL****LOT SPECIFIC VALUE ASSIGNMENT CARD****LEVEL 1****Haemoglobin A1c****UNITS****MEAN****RANGE**

NGSP ALIGNED

%

6.3

(5.3-7.3)

IFCC ALIGNED

mmol/mol

45

(38-52)

g/dL

0.61

(0.51-0.71)

mmol/L

0.38

(0.32-0.44)

**Total Haemoglobin**

g/dL

13.55

(11.17-15.93)

mmol/L

8.40

(6.93-9.87)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

**BECKMAN COULTER DxC CHEMISTRY SYSTEMS**

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4188



IVD

EC REP

Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New ZealandEMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands**USE BY 2021-08-31**

Issue 09. Revision date 09/17

# extendSURE™

**CONTROL****LOT SPECIFIC VALUE ASSIGNMENT CARD****LEVEL 2****Haemoglobin A1c****UNITS****MEAN****RANGE**

NGSP ALIGNED

%

10.7

(8.9-12.5)

IFCC ALIGNED

mmol/mol

93

(77-109)

g/dL

1.26

(1.04-1.48)

mmol/L

0.78

(0.64-0.92)

**Total Haemoglobin**

g/dL

13.41

(10.99-15.83)

mmol/L

8.32

(6.82-9.82)

MANUFACTURER'S INSTRUMENT/HAEMOGLOBIN A1c REAGENT ASSAY

**BECKMAN COULTER DxC CHEMISTRY SYSTEMS**

For use with Beckman Coulter HbA1c Reagent Kit Part Numbers A87905 and B36415

LOT

4188



IVD

EC REP

Canterbury Scientific Ltd  
71 Whiteleigh Avenue  
Addington, 8011  
Christchurch  
New ZealandEMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands**USE BY 2021-08-31**

Issue 09. Revision date 09/17

Printed to either side of  
225g card